


General	Good	-----
Nurse	Jannet Marigine	Yesterday: Anne Bargin
Activities	Walking- 20 min	30 Min
Sleep	7.5 Hours	Normal
Pressure	120/68	Less than 120/80 mmHG
Glucose	125 mg/dl	100 to 125 mg/dl
Oxygen	95%	95% to 100%
Medication	Blach 25mg, Blah 39mg, blah .5 mg	(Yesterday medication)

Call me
-Emergency-





Backend
Login

ID

Password

Log IN

Main page, main user

Patients

Nurses

Users

Data types



Backend

Main page, medical officer user

Patients

Nurses

Data types

Patients

	ID#	Data Sheet		Days in	Lines(month)	Pin Code -		Email
<input checked="" type="checkbox"/>	P056269	Download	View	57	5482	Show	Reset	Blah
<input checked="" type="checkbox"/>	P784664	Download	View	43	4531	Show	Reset	Blah
<input checked="" type="checkbox"/>	P668451	Download	View	28	2236	Show	Reset	Blah
New	-	Download	View	0	0	Show	Reset	Blah

ID number is entered

Add Patient

Archive

Sure, twice,
require note and
a check box.

Nurses

	ID#	First Name	Last name	Lines entered (Month)	Data Sheet (All what was entered by them-Month)		Password	-	Email
<input checked="" type="checkbox"/>	N2576	Blah	Blah	5546	Download	View	Show	Reset	Blah
<input checked="" type="checkbox"/>	N8875	Blah	Blah	6542	Download	View	Show	Reset	Blah
<input checked="" type="checkbox"/>	N5535	Blah	Blah	354654	Download	View	Show	Reset	Blah

ID number is entered

Add Nurse

Archive

Backend Users

	ID#	First Name	Last name	Activity		Password	-	Email
<input checked="" type="checkbox"/>	N2576	Blah	Blah	Download	View	Show	Reset	Blah
<input checked="" type="checkbox"/>	N8875	Blah	Blah	Download	View	Show	Reset	Blah
<input checked="" type="checkbox"/>	N5535	Blah	Blah	Download	View	Show	Reset	Blah

ID number is entered

Add Nurse

Archive

Data types

Move up and down, which will control the order it will be presented on the patient side device.

- If: answer is 0 to 99 then show "Good" in (green color)
- If: answer is 100 to 200 then show "Med" in (yellow color)

	ID#	Name	Type	Reference rules	Show in bar	Color
⬇️ <input checked="" type="checkbox"/>	DT257	Blah			Check Box	Hex Code
⬇️ <input checked="" type="checkbox"/>	DT875	Blah			Check Box	Hex Code
⬇️ <input checked="" type="checkbox"/>	DT535	Blah			Check Box	Hex Code

▼ Drop Down

Types:
numerical
Range
Text


Add

- Eat(numerical)
- Activities(text)
- Sleep(range)
- Nap(range)
- Shower(numerical)
- Move(numerical)
- Medication(text)
- Nurse name(text), reference

- would be yesterday's nurse)
- Pressure(numerical , require units)
- Oxygen(numerical , require units)
- Glycose(numerical , require units)

Delete

Nurse phone app



Patient Report Sheet



“Patient” Eyes Only

ID:
07-05-2021(Date)

(01-05-2021 00:00)

04-05-2021 00:00 Nurse ID: Nurse Name: Note ID
Type: Blah
Value:
Reference:

04-05-2021 00:00 Nurse ID: Nurse Name: Note ID
...

(06-05-2021 00:00)

END

End Date/Time, Unique Number for every print, Including used ID

Main

- The application is based but not limited to the presented design.
- The application, database and patients reports must be highly secure.
- The final system be reviewed by a third party to insure the quality of the security.
- Records can't be deleted.
- Records are sent automatically(end of month) by email to the Musandah email, and patient email. Including all the nurses one
- Some vitals signs on the patient side display will be removed after a certain amount of time(Assigned per type through backend).
- Two types of backend users, admin, and medical officer(who can't edit backend users).